


Novel Buoyancy Based Cell Selection: X-BACS™ Technology

Content originally prepared for ThermoGenesis by Dalip Sethi¹, Jon Ellis², Zelenia Contreras², John Perea², Stephen Truong², Jillian Miller¹
¹Cesca Therapeutics, Inc., Rancho Cordova, CA USA; ²ThermoGenesis Corp., Rancho Cordova, CA USA

Introduction


New developments in cell biology and genetic engineering have revolutionized the field of medicine with the greatest impact observed in the field of Immuno-oncology. Efficient isolation of desired cell populations is the cornerstone in the development of manufacturing and quality control processes for emerging cell therapies. We have developed a buoyancy-based cell selection method that isolates desired cell populations from a mixture of cell types, such as mononuclear cell (MNC) fraction preparations from whole blood. Using this method, we have efficiently recovered highly pure populations of T-cells (95%) with good yield (85%). The method is simple to execute with standard laboratory equipment.

Principle of Operation



Materials and Methods


Peripheral Blood (PB) units were purchased from BloodSource. MNC fractions were prepared using both Ficoll-Hypaque density gradient centrifugation or, as an alternative, the X-LAB® System.



Both MNC fractions were used to select CD3⁺ T-Cells using BACS technology. MNC fractions were incubated with anti-CD3 antibody for 30 minutes, followed by microbubble reagent for 20 minutes. Following incubation, the target and non-target fractions were separated using centrifugation.

Results

Whole blood units (n=4) were split among three (3) independent users, and each recovered MNCs using both Ficoll-Hypaque gradient centrifugation, as well as with the X-LAB System.



	X-LAB System			Ficoll Method		
	MNC (%) Recovery	PLT (%) Depletion	RBC (%) Depletion	MNC (%) Recovery	PLT (%) Depletion	RBC (%) Depletion
Mean	95.4	82.4	99.5	77.8	86.2	99.9
Std. Deviation	4.0	3.1	0.1	9.5	4.4	0.1
Std. Error of Mean	1.1	0.9	0.04	1.6	0.7	0.01
Lower 95% CI of mean	92.8	80.5	99.4	74.6	84.7	99.8
Upper 95% CI of mean	97.9	84.4	99.6	81.0	87.7	99.9

Figure 1. Comparative cell recoveries and depletions obtained using X-LAB Technology and Ficoll method.

The X-LAB system is a US class I medical device and medical device accessories. Should you decide that you want to use the product for any intended use other than that made by Corning then you accept that Corning does not make or support that claim and that you as the user are responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of your intended application. The X-BACS system is for research use only.
 ©2021 Corning Incorporated. All rights reserved.

CD3⁺ T-cells were selected from both MNC fractions using the X-BACS System. A mean CD3⁺ recovery of 83.5%, 96.3% viability, and 95.1% purity was obtained from Ficoll MNC fractions, while a mean CD3⁺ recovery of 89.0%, 98.3% viability, and 97.3% purity was obtained from X-LAB MNC fractions.




Figure 2. Comparative CD3⁺ selection data obtained from MNC fractions prepared with X-LAB Technology and Ficoll method.

The overall percent (%) CD3⁺ cell recovery was calculated by multiplying percent (%) MNC recovery with percent (%) CD3⁺ BACS recovery. In summary, better CD3⁺ selection was obtained using the X-LAB System followed by BACS compared to Ficoll method followed by BACS.




Figure 3. Comparative percent (%) CD3⁺ recovery data from starting material, whole blood, to purified CD3⁺ T-cells.

Conclusions

We have developed a buoyancy-based cell selection method that provides an efficient means of isolating a highly pure population of CD3⁺ T-cells (>95%) with good yield (>85%) and excellent viability (>95%).